

兔抗 CCL11 多克隆抗体

中文名称：兔抗 CCL11 多克隆抗体

英文名称： Anti-CCL11 rabbit polyclonal antibody

别名： C-C motif chemokine ligand 11 SCYA11

相关类别： 一抗

抗原： CCL11

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This antimicrobial gene is one of several chemokine genes clustered on the q-arm of chromosome 17. Chemokines form a superfamily of secreted proteins involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of the N-terminal cysteine residues of the mature peptide. This chemokine, a member of the CC subfamily, displays chemotactic activity for eosinophils, but not mononuclear cells or neutrophils. This eosinophil-specific chemokine is thought to be involved in eosinophilic inflammatory diseases such as atopic derma

	titis, allergic rhinitis, asthma and parasitic infections.
Applications:	ELISA, IHC
Name of antibody:	CCL11
Immunogen:	Synthetic peptide of human CCL11
Full name:	C-C motif chemokine ligand 11
Synonyms:	SCYA11
SwissProt:	P51671
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer;
IHC Recommend dilution:	20-100

