

## 兔抗 CCL25 多克隆抗体

中文名称：兔抗 CCL25 多克隆抗体

英文名称： Anti-CCL25 rabbit polyclonal antibody

别名： TECK; Ckb15; SCYA25

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CCL25

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This antimicrobial gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene displays chemotactic activity for dendritic cells, thymocytes, and activated macrophages but is inactive on peripheral blood lymphocytes and neutrophils. The product of this gene binds to chemokine receptor CCR9. Alternative splicing results in multiple t

	ranscript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CCL25
<b>Immunogen:</b>	Fusion protein of human CCL25
<b>Full name:</b>	chemokine (C-C motif) ligand 25
<b>Synonyms :</b>	TECK; Ckb15; SCYA25
<b>SwissProt:</b>	O15444
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human breast cancer
<b>IHC Recommend dilution:</b>	25-100

