

兔抗 DNAJC7 多克隆抗体

中文名称：兔抗 DNAJC7 多克隆抗体

英文名称： Anti-DNAJC7 rabbit polyclonal antibody

别名： DJ11, DJC7, TPR2, TTC2

抗原： DNAJC7

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the DNAJ heat shock protein 40 family of proteins that is characterized by two N-terminal tetratricopeptide repeat domains and a C-terminal DNAJ domain. This protein binds the chaperone proteins heat shock proteins 70 and 90 in an ATP-dependent manner and may function as a co-chaperone. Pseudogenes of this gene are found on chromosomes 1 and 6. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	DNAJC7
Immunogen:	Fusion protein of human DNAJC7

Full name:	DnaJ (Hsp40) homolog, subfamily C, member 7
Synonyms :	DJ11, DJC7, TPR2, TTC2
SwissProt:	Q99615
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	56 kDa
WB Positive control:	Human seminoma tissue
WB Recommended dilution:	500-2000



