

兔抗 CCN3 多克隆抗体

中文名称：兔抗 CCN3 多克隆抗体

英文名称： Anti-CCN3 rabbit polyclonal antibody

别名： cellular communication network factor 3; NOV; NOVh

相关类别： 一抗

抗原： CCN3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CNN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development. [provided by RefSeq, Feb 2009]
Applications:	ELISA, WB, IHC
Name of antibody:	CCN3
Immunogen:	Fusion protein of human CCN3
Full name:	cellular communication network factor 3
Synonyms:	NOV; NOVh; IBP-9; IGFBP9; IGFBP-9

SwissProt:	P48745
IHC positive control:	Human ovarian cancer and Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	39 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	500-2000



