

兔抗 CCND1(Phospho-Ser90) 多克隆抗体

中文名称：兔抗 CCND1(Phospho-Ser90) 多克隆抗体

英文名称： Anti-CCND1(Phospho-Ser90) rabbit polyclonal antibody

别名： BCL1; PRAD1; U21B31; D11S287E

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CCND1(Phospho-Ser90)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance throughout the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with and functions as a regulatory subunit of CDK4 or CDK6, whose activity is required for cell cycle G1/S transition.
Applications:	WB

Name of antibody:	CCND1(Phospho-Ser90)
Immunogen:	Peptide sequence around phosphorylation site of serine 90(F-L-S(p)-L-E) derived from Human Cyclin D1 .
Full name:	cyclin D1
Synonyms :	BCL1; PRAD1; U21B31; D11S287E
SwissProt:	P24385
WB Predicted band size:	34 kDa
WB Positive control:	K562 cells lysates
WB Recommended dilution:	500-1000

