

兔抗 DND1 多克隆抗体

中文名称：兔抗 DND1 多克隆抗体

英文名称：Anti-DND1 rabbit polyclonal antibody

别名：RBMS4

相关类别：一抗

抗原：DND1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that binds to microRNA-targeting sequences of mRNAs, inhibiting microRNA-mediated repression. Reduced expression of this gene has been implicated in tongue squamous cell carcinoma. Two pseudogenes of this gene are located on the long arm of chromosome 17.
Applications:	ELISA, IHC
Name of antibody:	DND1
Immunogen:	Synthetic peptide of human DND1
Full name:	DND microRNA-mediated repression inhibitor 1
Synonyms:	RBMS4

SwissProt:	Q8IYX4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;Human breast cancer
IHC Recommend dilution:	20-100



