

兔抗 DOK5 多克隆抗体

中文名称：兔抗 DOK5 多克隆抗体

英文名称： Anti-DOK5 rabbit polyclonal antibody

别名： C20orf180

储存： 冷冻（-20℃）

抗原： DOK5

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the DOK family of membrane proteins, which are adapter proteins involved in signal transduction. The encoded protein interacts with phosphorylated receptor tyrosine kinases to mediate neurite outgrowth and activation of the MAP kinase pathway. In contrast to other DOK family proteins, this protein does not interact with RASGAP.
Applications:	ELISA, IHC
Name of antibody:	DOK5
Immunogen:	Fusion protein of human DOK5
Full name:	docking protein 5

Synonyms :	C20orf180
SwissProt:	Q9P104
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	10-50

