

兔抗 CD109 多克隆抗体

中文名称：兔抗 CD109 多克隆抗体

英文名称： Anti-CD109 rabbit polyclonal antibody

别名： CPAMD7; RP11-525G3.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CD109

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the alpha2-macroglobulin/complement superfamily. The encoded GPI-linked glycoprotein is found on the cell surface of platelets, activated T-cells, and endothelial cells. The protein binds to and negatively regulates signaling of transforming growth factor beta (TGF-beta). Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CD109
Immunogen:	Synthetic peptide of human CD109

Full name:	CD109 molecule
Synonyms :	CPAMD7; RP11-525G3.1
SwissProt:	Q6YHK3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

