

兔抗 DSTYK 多克隆抗体

中文名称：兔抗 DSTYK 多克隆抗体

英文名称： Anti-DSTYK rabbit polyclonal antibody

别名： RIP5; RIPK5; CAKUT1; DustyPK; HDCMD38P

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DSTYK

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a dual serine/threonine and tyrosine protein kinase which is expressed in multiple tissues. It is thought to function as a regulator of cell death. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	DSTYK
Immunogen:	Fusion protein of human DSTYK
Full name:	dual serine/threonine and tyrosine protein kinase
Synonyms:	RIP5; RIPK5; CAKUT1; DustyPK; HDCMD38P
SwissProt:	Q6XUX3

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	105 kDa
WB Positive control:	823 and hepg2 cells
WB Recommended dilution:	200-1000



