

## ASMTL 抗原（重组蛋白）

中文名称： ASMTL 抗原（重组蛋白）

英文名称： ASMTL Antigen (Recombinant Protein)

别名： acetylserotonin O-methyltransferase like; ASTML; ASMTLX; ASMTLY

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 422-621 amino acids of human ASMTL

### 技术规格

<b>Full name:</b>	acetylserotonin O-methyltransferase like
<b>Synonyms:</b>	ASTML; ASMTLX; ASMTLY
<b>Swissprot:</b>	O95671
<b>Gene Accession:</b>	BC002508
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentation (maf) protein of <i>Bacillus subtilis</i> and to orfE of <i>Escherichia coli</i> , while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene.