

ATAD3A 抗原（重组蛋白）

中文名称： ATAD3A 抗原（重组蛋白）

英文名称： ATAD3A Antigen (Recombinant Protein)

别名： ATPase family, AAA domain containing 3A

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 435-634 amino acids of human ATAD3A

技术规格

Full name:	ATPase family, AAA domain containing 3A
Swissprot:	Q9NVI7
Gene Accession:	BC007803
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	ATAD3A and ATAD3B (MIM 612317) are mitochondrial membrane proteins that contribute to the stabilization of large mitochondrial DNA (mtDNA)-protein complexes called nucleoids (He et al., 2007 [PubMed 17210950]).

mlbio 梅联生物
Good elisakit producers