

兔抗 CDADC1 多克隆抗体

中文名称：兔抗 CDADC1 多克隆抗体

英文名称： Anti-CDADC1 rabbit polyclonal antibody

别名： cytidine and dCMP deaminase domain containing 1; NYD-SP15; bA103J18.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CDADC1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

NYD-SP15 (testis development protein NYD-SP15), also known as CDADC1 (cytidine and dCMP deaminase domain containing 1), is a widely expressed protein with predominant expression in the testis, liver, spleen, kidney, thymus and placenta. NYD-SP15 is 514 amino acids in length and belongs to the cytidine and deoxycytidylate deaminase family. It is developmentally regulated with higher expression in adult testis than fetal testis and is believed to participate in spermatogenesis and testicular development. This suggests that NYD-SP15 may be a determining factor in male infertility. Due to alternative

	splicing events, four transcript variants exist for NYD-S P15.
Applications:	ELISA, IHC
Name of antibody:	CDADC1
Immunogen:	Full length fusion protein
Full name:	cytidine and dCMP deaminase domain containing 1
Synonyms:	NYD-SP15; bA103J18.1
SwissProt:	Q9BWV3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	30-150

