

ATF7 抗原（重组蛋白）

中文名称： ATF7 抗原（重组蛋白）

英文名称： ATF7 Antigen (Recombinant Protein)

别名： ATFA

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	activating transcription factor 7
Synonyms:	ATFA
Swissprot:	P17544
Gene Accession:	BC042363
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Plays important functions in early cell signaling. Binds the cAMP response element (CRE) (consensus: 5'-GTGACGT[AG][AG]-3'), a sequence present in many viral and cellular promoters. Activator of the NF-ELAM1/delta-A site of the E-selectin promoter. Has no intrinsic transcriptional activity, but activates transcription on formation of JUN or FOS heterodimers. Also can bind TRE promoter sequences when hetero

dimerized with members of the JUN family.