

兔抗 CDC25A(Phospho-Ser178) 多克隆抗体

中文名称：兔抗 CDC25A(Phospho-Ser178) 多克隆抗体

英文名称： Anti-CDC25A(Phospho-Ser178) rabbit polyclonal antibody

别名： CDC25A2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CDC25A(Phospho-Ser178)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CDC25A is a member of the CDC25 family of phosphatases. CDC25A is required for progression from G1 to the S phase of the cell cycle. It activates the cyclin-dependent kinase CDC2 by removing two phosphate groups. CDC25A is specifically degraded in response to DNA damage, which prevents cells with chromosomal abnormalities from progressing through cell division.
Applications:	WB
Name of antibody:	CDC25A(Phospho-Ser178)
Immunogen:	Peptide sequence around phosphorylation site of s

	erine 178(Q-N-S(p)-A-P) derived from Human CDC2 5A .
Full name:	cell division cycle 25A
Synonyms :	CDC25A2
SwissProt:	P30304
WB Predicted band size:	59 kDa
WB Positive control:	A549 cells lysates treated with UV
WB Recommended dilution:	500-1000

