

兔抗 CDC37L1 多克隆抗体

中文名称：兔抗 CDC37L1 多克隆抗体

英文名称： Anti-CDC37L1 rabbit polyclonal antibody

别名： HARC; CDC37B; RP11-6J24.5

相关类别： 一抗

抗原： CDC37L1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding
Applications:	ELISA, IHC
Name of antibody:	CDC37L1
Immunogen:	Fusion protein of human CDC37L1
Full name:	cell division cycle 37-like 1
Synonyms :	HARC; CDC37B; RP11-6J24.5
SwissProt:	Q7L3B6
ELISA Recommended dilution:	2000-10000

IHC positive control:	Human cervical cancer and human colon cancer
IHC Recommend dilution:	100-300

