

## 兔抗 CDC42BPB 多克隆抗体

中文名称：兔抗 CDC42BPB 多克隆抗体

英文名称： Anti-CDC42BPB rabbit polyclonal antibody

别名： MRCKB

抗原： CDC42BPB

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the serine/threonine protein kinase family. The encoded protein contains a Cdc42/Rac-binding p21 binding domain resembling that of PAK kinase. The kinase domain of this protein is most closely related to that of myotonic dystrophy kinase-related ROK. Studies of the similar gene in rat suggested that this kinase may act as a downstream effector of Cdc42 in cytoskeletal reorganization.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CDC42BPB
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MRC

	KB.
<b>Full name:</b>	CDC42 binding protein kinase beta (DMPK-like)
<b>Synonyms :</b>	MRCKB
<b>SwissProt:</b>	Q9Y5S2
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	194 kDa
<b>WB Positive control:</b>	COLO cells lysate
<b>WB Recommended dilution:</b>	500-3000



