

兔抗 CDCA5 多克隆抗体

中文名称: 兔抗 CDCA5 多克隆抗体

英文名称: Anti-CDCA5 rabbit polyclonal antibody

别名: SORORIN

抗原: CDCA5

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Regulator of sister chromatid cohesion in mitosis stabilizing cohesin complex association with chromatin. May antagonize the action of WAPL which stimulates cohesin dissociation from chromatin. Cohesion ensures that chromosome partitioning is accurate in both meiotic and mitotic cells and plays an important role in DNA repair. Required for efficient DNA double-stranded break repair.
Applications:	ELISA, IHC
Name of antibody:	CDCA5
Immunogen:	Synthetic peptide of human CDCA5

Full name:	cell division cycle associated 5
Synonyms:	SORORIN
SwissProt:	Q96FF9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and human thyroid cancer
IHC Recommend dilution:	25-100



