

## 兔抗 DVL1 多克隆抗体

中文名称: 兔抗 DVL1 多克隆抗体

英文名称: Anti-DVL1 rabbit polyclonal antibody

别名: DVL; DRS2; DVL1L1; DVL1P1

抗原: DVL1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

**Background:**

DVL1, the human homolog of the Drosophila dishevelled gene (dsh) encodes a cytoplasmic phosphoprotein that regulates cell proliferation, acting as a transducer molecule for developmental processes, including segmentation and neuroblast specification. DVL1 is a candidate gene for neuroblastomatous transformation. The Schwartz-Jampel syndrome and Charcot-Marie-Tooth disease type 2A have been mapped to the same region as DVL1. The phenotypes of these diseases may be consistent with defects which might be expected from aberrant expression of a DVL gene during development.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DVL1
<b>Immunogen:</b>	Synthetic peptide of human DVL1
<b>Full name:</b>	dishevelled segment polarity protein 1
<b>Synonyms :</b>	DVL; DRS2; DVL1L1; DVL1P1
<b>SwissProt:</b>	O14640
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100



