

兔抗 DVL3 多克隆抗体

中文名称：兔抗 DVL3 多克隆抗体

英文名称： Anti-DVL3 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： DVL3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.
Applications:	ELISA, IHC
Name of antibody:	DVL3
Immunogen:	Synthetic peptide of human DVL3
Full name:	dishevelled segment polarity protein 3
SwissProt:	Q92997
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and Human brain
IHC Recommend dilution:	25-100

