

## 兔抗 DYDC2 多克隆抗体

中文名称：兔抗 DYDC2 多克隆抗体

英文名称：Anti-DYDC2 rabbit polyclonal antibody

相关类别：一抗

抗 原：DYDC2

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of proteins that contains a DPY30 domain. This gene locus overlaps with a closely related gene on the opposite strand. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DYDC2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	DPY30 domain containing 2
<b>SwissProt:</b>	Q96IM9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

**mlbio** 梅联生物  
Good elisakit producers

