

BARX1 抗原（重组蛋白）

中文名称： BARX1 抗原（重组蛋白）

英文名称： BARX1 Antigen (Recombinant Protein)

别名： BARX homeobox 1

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 155-254 amino acids of human BARX1

技术规格

Full name:	BARX homeobox 1
Swissprot:	Q9HBU1
Gene Accession:	BC009458
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the Bar subclass of homeobox transcription factors. Studies of the mouse and chick homolog suggest the encoded protein may play a role in developing teeth and craniofacial mesenchyme of neural crest origin. The protein may also be associated with differentiation of stomach epithelia.