

兔抗 DYNAP 多克隆抗体

中文名称：兔抗 DYNAP 多克隆抗体

英文名称：Anti-DYNAP rabbit polyclonal antibody

别名：C18orf26

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DYNAP

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Plays a role in the regulation of cell proliferation. Promotes activation of the AKT1 signaling pathway. Promotes phosphorylation of AKT1 at 'Ser-473'.
Applications:	ELISA, IHC
Name of antibody:	DYNAP
Immunogen:	Synthetic peptide of human DYNAP
Full name:	dynactin associated protein
Synonyms:	C18orf26
SwissProt:	Q8N1N2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain

IHC Recommend dilution:

30-150

