

兔抗 DYNLRB1 多克隆抗体

中文名称：兔抗 DYNLRB1 多克隆抗体

英文名称：Anti-DYNLRB1 rabbit polyclonal antibody

别名：dynein light chain roadblock-type 1; BLP; BITH; DNCL2A; DNLC2A; ROBLD1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DYNLRB1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene is a member of the roadblock dynein light chain family. The encoded cytoplasmic protein is capable of binding intermediate chain proteins, interacts with transforming growth factor-beta, and has been implicated in the regulation of actin modulating proteins. Upregulation of this gene has been associated with hepatocellular carcinomas, suggesting that this gene may be involved in tumor progression. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 12 and 18. [provided by RefSeq, Aug 2013]

Applications:	ELISA, IHC
Name of antibody:	DYNLRB1
Immunogen:	Fusion protein of human DYNLRB1
Full name:	dynein light chain roadblock-type 1
Synonyms:	BLP; BITH; DNCL2A; DNLC2A; ROBLD1
SwissProt:	Q9NP97
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and Human tonsil
IHC Recommend dilution:	150-300



