

兔抗 DYRK1B 多克隆抗体

中文名称：兔抗 DYRK1B 多克隆抗体

英文名称： Anti-DYRK1B rabbit polyclonal antibody

别名： MIRK; AOMS3

抗原： DYRK1B

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of a family of nuclear-localized protein kinases. The encoded protein participates in the regulation of the cell cycle. Expression of this gene may be altered in tumor cells, and mutations in this gene were found to cause abdominal obesity-metabolic syndrome 3. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	DYRK1B
Immunogen:	Synthesized peptide derived from internal of human DYRK1B.

Full name:	dual-specificity tyrosine-(Y)-phosphorylation regulated kinase 1B
Synonyms :	MIRK; AOMS3
SwissProt:	Q9Y463
WB Predicted band size:	69 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

