

BCAN 抗原（重组蛋白）

中文名称： BCAN 抗原（重组蛋白）

英文名称： BCAN Antigen (Recombinant Protein)

别名： brevican; BEHAB; CSPG7

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 712-911 amino acids of human BCAN

技术规格

Full name:	brevican
Synonyms:	BEHAB; CSPG7
Swissprot:	Q96GW7
Gene Accession:	BC022938
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants.

mlbio 梅联生物
Good elisakit producers