

兔抗 E2F1(Phospho-Thr433) 多克隆抗体

中文名称：兔抗 E2F1(Phospho-Thr433) 多克隆抗体

英文名称： Anti-E2F1(Phospho-Thr433) rabbit polyclonal antibody

别名： RBP3; E2F-1; RBAP1; RBBP3

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： E2F1(Phospho-Thr433)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the E2F family of transcription factors. The E2F family plays a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of the transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserved domains found in most members of the family.
Applications:	WB
Name of antibody:	E2F1(Phospho-Thr433)
Immunogen:	Peptide sequence around phosphorylation site of Threonine 433 (D-L-T(p)-P-L) derived from Human E2F1.

Full name:	E2F transcription factor 1
Synonyms :	RBP3; E2F-1; RBAP1; RBBP3
SwissProt:	Q01094
WB Predicted band size:	47 kDa
WB Positive control:	Mouse brain tissue lysates
WB Recommended dilution:	500-1000

