

兔抗 E2F1 多克隆抗体

中文名称：兔抗 E2F1 多克隆抗体

英文名称：Anti-E2F1 rabbit polyclonal antibody

别名：RBP3; E2F-1; RBAP1; RBBP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：E2F1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the E2F family of transcription factors. The E2F family plays a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of the transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserved domains found in most members of the family. These domains include a DNA binding domain, a dimerization domain which determines interaction with the differentiation regulated transcription factor proteins (DP), a transactivation domain enriched in acidic amino acids, and a tumor suppressor

	essor protein association domain which is embedded within the transactivation domain.
Applications:	ELISA, WB
Name of antibody:	E2F1
Immunogen:	Synthetic peptide of human E2F1
Full name:	E2F transcription factor 1
Synonyms :	RBP3; E2F-1; RBAP1; RBBP3
SwissProt:	Q01094
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	47 kDa
WB Positive control:	Human fetal muscle tissue
WB Recommended dilution:	200-1000

