

兔抗 CDHR4 多克隆抗体

中文名称: 兔抗 CDHR4 多克隆抗体

英文名称: Anti-CDHR4 rabbit polyclonal antibody

别名: cadherin related family member 4; CDH29; PRO34300

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: CDHR4

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types (By similarity).
Applications:	ELISA, IHC
Name of antibody:	CDHR4
Immunogen:	Synthetic peptide of human CDHR4
Full name:	cadherin related family member 4
Synonyms:	CDH29; PRO34300
SwissProt:	A6H8M9

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

