

兔抗 E2F4 多克隆抗体

中文名称：兔抗 E2F4 多克隆抗体

英文名称：Anti-E2F4 rabbit polyclonal antibody

别名：E2F-4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：E2F4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the E2F family of transcription factors. The E2F family play s a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of t he transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserv ed domains found in most members of the family. Th ese domains include a DNA binding domain, a dimeri zation domain which determines interaction with the differentiation regulated transcription factor proteins (DP), a transactivation domain enriched in acidic amino

	acids, and a tumor suppressor protein association domain which is embedded within the transactivation domain. This protein binds to all three of the tumor suppressor proteins pRB, p107 and p130, but with higher affinity to the last two.
Applications:	ELISA, IHC
Name of antibody:	E2F4
Immunogen:	Synthetic peptide of human E2F4
Full name:	E2F transcription factor 4, p107/p130-binding
Synonyms :	E2F-4
SwissProt:	Q16254
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	50-200

