

兔抗 ECI1 多克隆抗体

中文名称：兔抗 ECI1 多克隆抗体

英文名称：Anti-ECI1 rabbit polyclonal antibody

别名：DCI

抗原：ECI1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the hydratase/isomerase superfamily. The protein encoded is a key mitochondrial enzyme involved in beta-oxidation of unsaturated fatty acids. It catalyzes the transformation of 3-cis and 3-trans-enoyl-CoA esters arising during the stepwise degradation of cis-, mono-, and polyunsaturated fatty acids to the 2-trans-enoyl-CoA intermediates. Alternatively spliced transcript variants have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	ECI1
Immunogen:	Fusion protein of human ECI1

Full name:	enoyl-CoA delta isomerase 1
Synonyms :	DCI
SwissProt:	P42126
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	33 kDa
WB Positive control:	Human fetal liver tissue , MCF-7 cell , 293T cell , Hela cell , A549 cell
WB Recommended dilution:	500-2000



