

兔抗 ECRG4 多克隆抗体

中文名称：兔抗 ECRG4 多克隆抗体

英文名称：Anti-ECRG4 rabbit polyclonal antibody

别名：ECRG4 augurin precursor; C2orf40

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ECRG4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Probable hormone that may attenuate cell proliferation and induce senescence of oligodendrocyte and neural precursor cells in the central nervous system (By similarity). ECRG4-induced senescence is characterized by G1 arrest, RB1 dephosphorylation and accelerated CCND1 and CCND3 proteasomal degradation (By similarity).
Applications:	ELISA, IHC
Name of antibody:	ECRG4
Immunogen:	Synthetic peptide of human ECRG4
Full name:	ECRG4 augurin precursor
Synonyms:	C2orf40

SwissProt:	Q9H1Z8
IHC positive control:	Human thyroid cancer and Human brain
IHC Recommend dilution:	25-100

