

## BLVRA 抗原（重组蛋白）

中文名称： BLVRA 抗原（重组蛋白）

英文名称： BLVRA Antigen (Recombinant Protein)

别名： BVR; BLVR; BVRA

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	biliverdin reductase A
<b>Synonyms:</b>	BVR; BLVR; BVRA
<b>Swissprot:</b>	P53004
<b>Gene Accession:</b>	BC008456
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene belongs to the biliverdin reductase family, members of which catalyze the conversion of biliverdin to bilirubin in the presence of NADPH or NADH. Mutations in this gene are associated with hyperbiliverdinemia. Alternatively spliced transcript variants have been found for this gene.