

兔抗 EDEM1 多克隆抗体

中文名称：兔抗 EDEM1 多克隆抗体

英文名称：Anti-EDEM1 rabbit polyclonal antibody

别名：EDEM

相关类别：一抗

抗原：EDEM1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Extracts misfolded glycoproteins, but not glycoproteins undergoing productive folding, from the calnexin cycle. It is directly involved in endoplasmic reticulum-associated degradation (ERAD) and targets misfolded glycoproteins for degradation in an N-glycan-independent manner, probably by forming a complex with SEL1L. It has low mannosidase activity, catalyzing mannose trimming from Man8GlcNAc2 to Man7GlcNAc2.
Applications:	ELISA, IHC
Name of antibody:	EDEM1

Immunogen:	Synthetic peptide of human EDEM1
Full name:	ER degradation enhancing alpha-mannosidase like protein 1
Synonyms:	EDEM
SwissProt:	Q92611
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;Human gastric cancer
IHC Recommend dilution:	10-50



