

## 兔抗 EDEM2 多克隆抗体

中文名称: 兔抗 EDEM2 多克隆抗体

英文名称: Anti-EDEM2 rabbit polyclonal antibody

别名: C20orf31; C20orf49; bA4204.1

储存: 冷冻 (-20℃)

抗原: EDEM2

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	In the endoplasmic reticulum (ER), misfolded proteins are retrotranslocated to the cytosol and degraded by the proteasome in a process known as ER-associated degradation (ERAD). EDEM2 belongs to a family of proteins involved in ERAD of glycoproteins (Mast et al., 2005 [PubMed 15537790]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EDEM2
<b>Immunogen:</b>	Fusion protein of human EDEM2
<b>Full name:</b>	ER degradation enhancing alpha-mannosidase like protein 2

<b>Synonyms:</b>	C20orf31; C20orf49; bA4204.1
<b>SwissProt:</b>	Q9BV94
<b>ELISA Recommended dilution:</b>	500-1000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

