

兔抗 EDEM3 多克隆抗体

中文名称：兔抗 EDEM3 多克隆抗体

英文名称：Anti-EDEM3 rabbit polyclonal antibody

别名：ER degradation enhancing alpha-mannosidase like protein 3; C1orf22

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EDEM3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Quality control in the endoplasmic reticulum (ER) ensures that only properly folded proteins are retained in the cell through recognition and degradation of misfolded or un assembled proteins. EDEM3 belongs to a group of proteins that accelerate degradation of misfolded glycoproteins in the ER (Hirao et al., 2006 [PubMed 16431915]).
Applications:	ELISA, IHC
Name of antibody:	EDEM3
Immunogen:	Synthetic peptide of human EDEM3
Full name:	ER degradation enhancing alpha-mannosidase like protein 3

Synonyms:	C1orf22
SwissProt:	Q9BZQ6
ELISA Recommended dilution:	500-1000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

