

## BRMS1L 抗原（重组蛋白）

中文名称： BRMS1L 抗原（重组蛋白）

英文名称： BRMS1L Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

别 名： breast cancer metastasis-suppressor 1-like; BRMS1

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 224-323 amino acids of human BRMS1L

### 技术规格

<b>Full name:</b>	breast cancer metastasis-suppressor 1-like
<b>Synonyms:</b>	BRMS1
<b>Swissprot:</b>	Q5PSV4
<b>Gene Accession:</b>	BC006250
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene shows sequence similarity to the human breast carcinoma metastasis suppressor (BRMS1) protein and the mammalian Sds3 (suppressor of defective silencing 3) proteins. This protein is a component of the mSin3a family of histone deacetylase complexes (HDAC).