

## 兔抗 EDN1 多克隆抗体

中文名称：兔抗 EDN1 多克隆抗体

英文名称： Anti-EDN1 rabbit polyclonal antibody

别名： ET1; PPET1; HDLCQ7

相关类别： 一抗

抗原： EDN1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is proteolytically processed to release a secreted peptide termed endothelin 1. This peptide is a potent vasoconstrictor and is produced by vascular endothelial cells. Endothelin 1 also can affect the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EDN1
<b>Immunogen:</b>	Fusion protein of human EDN1
<b>Full name:</b>	endothelin 1
<b>Synonyms :</b>	ET1; PPET1; HDLCQ7

<b>SwissProt:</b>	P05305
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

