

BRPF3 抗原（重组蛋白）

中文名称： BRPF3 抗原（重组蛋白）

英文名称： BRPF3 Antigen (Recombinant Protein)

别名： bromodomain and PHD finger containing 3

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 1006-1205 amino acids of human BRPF3

技术规格

Full name:	bromodomain and PHD finger containing 3
Swissprot:	Q9ULD4
Gene Accession:	BC117387
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity.