

## 兔抗 EDN2 多克隆抗体

中文名称：兔抗 EDN2 多克隆抗体

英文名称： Anti-EDN2 rabbit polyclonal antibody

别名： ET2; PPET2

相关类别： 一抗

抗原： EDN2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the endothelin protein family of secretory vasoconstrictive peptides. The preproprotein is processed to a short mature form which functions as a ligand for the endothelin receptors that initiate intracellular signaling events. This gene product is involved in a wide range of biological processes, such as hypertension and ovulation.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	EDN2
<b>Immunogen:</b>	Fusion protein of human EDN2
<b>Full name:</b>	endothelin 2

<b>Synonyms :</b>	ET2; PPET2
<b>SwissProt:</b>	P20800
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and human colon cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	20 kDa
<b>WB Positive control:</b>	Human normal right lung and normal left lung tissue
<b>WB Recommended dilution:</b>	500-2000



