

BSA 抗原

中文名称: BSA 抗原

英文名称: BSA Antigen

相关类别: 抗原

储 存: 冷冻 (-20℃) 避光

概述

Full length protein of BSA

技术规格

Full name:	Bovine serum albumin
Swissprot:	P02769
Gene Accession:	NP_851335
Background:	Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca ²⁺ , Na ⁺ , K ⁺ , fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.