

兔抗 CDKAL1 多克隆抗体

- 中文名称：兔抗 CDKAL1 多克隆抗体
- 英文名称：Anti-CDKAL1 rabbit polyclonal antibody
- 别名：CDK5 regulatory subunit associated protein 1 like 1
- 相关类别：一抗
- 抗原：CDKAL1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.
Applications:	ELISA, IHC
Name of antibody:	CDKAL1
Immunogen:	Full length fusion protein
Full name:	CDK5 regulatory subunit associated protein 1 like 1
SwissProt:	Q5VV42
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300

