

## 兔抗 EEF1G 多克隆抗体

中文名称：兔抗 EEF1G 多克隆抗体

英文名称： Anti-EEF1G rabbit polyclonal antibody

别名： EF1G, GIG35

抗原： EEF1G

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains a n N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases. Probably plays a role in anchoring the complex to other cellular components.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	EEF1G
<b>Immunogen:</b>	Fusion protein of human EEF1G

<b>Full name:</b>	eukaryotic translation elongation factor 1 gamma
<b>Synonyms :</b>	EF1G, GIG35
<b>SwissProt:</b>	P26641
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Human liver cancer tissue
<b>WB Recommended dilution:</b>	500-2000



