

## 兔抗 EFEMP1 多克隆抗体

中文名称：兔抗 EFEMP1 多克隆抗体

英文名称： Anti-EFEMP1 rabbit polyclonal antibody

别名： DHRD; DRAD; FBNL; MLVT; MTLV; S1-5; FBLN3; FIBL-3

抗原： EFEMP1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a member of the fibulin family of extracellular matrix glycoproteins. Like all members of this family, the encoded protein contains tandemly repeated epidermal growth factor-like repeats followed by a C-terminus fibulin-type domain. This gene is upregulated in malignant gliomas and may play a role in the aggressive nature of these tumors. Mutations in this gene are associated with Doyme honeycomb retinal dystrophy. Alternatively spliced transcript variants that encode the same protein have been described.

<b>Applications:</b>	WB
<b>Name of antibody:</b>	EFEMP1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human EFEMP1.
<b>Full name:</b>	EGF containing fibulin-like extracellular matrix protein 1
<b>Synonyms :</b>	DHRD; DRAD; FBNL; MLVT; MTLV; S1-5; FBLN3; FBL-3
<b>SwissProt:</b>	Q12805
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Rat lung tissue lysate
<b>WB Recommended dilution:</b>	500-3000

