

## 兔抗 EFEMP2 多克隆抗体

中文名称：兔抗 EFEMP2 多克隆抗体

英文名称： Anti-EFEMP2 rabbit polyclonal antibody

别名： MBP1; UPH1; FBLN4; ARCL1B

抗原： EFEMP2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

A large number of extracellular matrix proteins have been found to contain variations of the epidermal growth factor (EGF) domain and have been implicated in functions as diverse as blood coagulation, activation of complement and determination of cell fate during development. The protein encoded by this gene contains four EGF2 domains and six calcium-binding EGF2 domains. This gene is necessary for elastic fiber formation and connective tissue development. Defects in this gene are cause of an autosomal recessive cutis laxa syndrome.

<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	EFEMP2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human EFEMP2.
<b>Full name:</b>	EGF containing fibulin-like extracellular matrix protein 2
<b>Synonyms :</b>	MBP1; UPH1; FBLN4; ARCL1B
<b>SwissProt:</b>	O95967
<b>IHC positive control:</b>	Human placenta tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	49 kDa
<b>WB Positive control:</b>	K562 cells and NIH/3T3 cells lysates
<b>WB Recommended dilution:</b>	500-3000
<b>IF positive control:</b>	A549 cells
<b>IF Recommend dilution:</b>	100-500



