

兔抗 CDO1 多克隆抗体

中文名称：兔抗 CDO1 多克隆抗体

英文名称： Anti-CDO1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： CDO1

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CDO1 (cysteine dioxygenase, type I) is a 200 amino acid protein that belongs to the cysteine dioxygenase family and is involved in organosulfur biosynthesis. Existing as a monomer and expressed at high levels in liver and placenta and at lower levels in brain, pancreas and heart, CDO1 functions as a dioxygenase that uses iron and zinc as cofactors to catalyze the conversion of L-cysteine and oxygen to 3-sulfinoalanine.
Applications:	ELISA, IHC
Name of antibody:	CDO1
Immunogen:	Fusion protein of human CDO1
Full name:	cysteine dioxygenase type 1

SwissProt:	Q16878
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100

