

CABP5 抗原（重组蛋白）

中文名称： CABP5 抗原（重组蛋白）

英文名称： CABP5 Antigen (Recombinant Protein)

别名： CABP3

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	calcium binding protein 5
Synonyms:	CABP3
Swissprot:	Q9NP86
Gene Accession:	BC126133
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The product of this gene belongs to a subfamily of calcium binding proteins, which share similarity to calmodulin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Expression of this gene is retina-specific. The mouse homolog of this protein has been shown to express