

CALML5 抗原（重组蛋白）

中文名称： CALML5 抗原（重组蛋白）

英文名称： CALML5 Antigen (Recombinant Protein)

别名： CLSP

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	calmodulin like 5
Synonyms:	CLSP
Swissprot:	Q9NZT1
Gene Accession:	BC039172
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a novel calcium binding protein expressed in the epidermis and related to the calmodulin family of calcium binding proteins. Functional studies with recombinant protein demonstrate it does bind calcium and undergoes a conformational change when it does so. Abundant expression is detected only in reconstructed epidermis and is restricted to differentiating keratinocytes. In addition, it can associate with transglutaminase.

minase 3, shown to be a key enzyme in the terminal differentiation of keratinocytes.