

## 兔抗 CDRT15L2 多克隆抗体

中文名称：兔抗 CDRT15L2 多克隆抗体

英文名称：Anti-CDRT15L2 rabbit polyclonal antibody

别名：CMT1A duplicated region transcript 15-like 2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CDRT15L2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CDRT15L2
<b>Immunogen:</b>	Synthetic peptide of human CDRT15L2
<b>Full name:</b>	CMT1A duplicated region transcript 15-like 2
<b>SwissProt:</b>	A8MXV6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

