

兔抗 CDX1 多克隆抗体

中文名称：兔抗 CDX1 多克隆抗体

英文名称： Anti-CDX1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： CDX1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity.
Applications:	ELISA, IHC
Name of antibody:	CDX1
Immunogen:	Fusion protein of human CDX1
Full name:	caudal type homeobox 1
SwissProt:	P47902
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

