

## 兔抗 EFNB1 多克隆抗体

中文名称：兔抗 EFNB1 多克隆抗体

英文名称：Anti-EFNB1 rabbit polyclonal antibody

级别：CFND; CFNS; EFB1; EFL3; EPLG2; Elk-L; LERK2

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：EFNB1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

相关类别：一抗

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a type I membrane protein and a ligand of Eph-related receptor tyrosine kinases. It may play a role in cell adhesion and function in the development or maintenance of the nervous system.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EFNB1
<b>Immunogen:</b>	Synthetic peptide of human EFNB1
<b>Full name:</b>	ephrin-B1
<b>Synonyms:</b>	CFND; CFNS; EFB1; EFL3; EPLG2; Elk-L; LERK2
<b>SwissProt:</b>	P98172

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	30-150

